

The difference makes it happen.

MiniMizer Extra Gastric Band

MiniMIZER
EXTRA





The difference makes it happen.
MiniMizer Extra Gastric Band



MiniMIZER®
EXTRA

Fixation without
plication.






MiniMIZER EXTRA

elastic loops to avoid band dislocation


The major differences:

- 1 direct fixation without plication:
reduced slippage
- 2 two closing positions:
less postoperative adjustments
- 3 inner compartment
integrated in band:
precise and concentric filling



first closing position 10.5 cm (\varnothing 31 mm)
second closing position 9.75 cm (\varnothing 26 mm)


precise and concentric filling



MiniMIZER EXTRA



- 4 pre-shaped gastric band for easy positioning
- 5 small port (ø 13 mm) for sub-xyphoidal placement
- 6 medium port (ø 18 mm) for sub-costal or rectal sheath placement



MiniMIZER EXTRA



with loops



 **HOSPI MEDICAL**
Bariatric Solutions

Hospimedical GmbH
Gewerbegebiet Egelsee
CH - 8259 Kaltenbach

Phone: +41 - 52 741 - 65 00

Fax: +41 - 52 741 - 65 02

info@hospimedical.com

www.hospimedical.com