

The difference makes it happen.

MiniMizer Extra Gastric Band

*Mini*MIZER
EXTRA





The difference makes it happen.
MiniMizer Extra Gastric Band



MiniMIZER®
EXTRA

Fixation without
plication.






MiniMIZER EXTRA

elastic loops to avoid band dislocation

The major differences:

- 1 direct fixation without plication:
reduced slippage
- 2 two closing positions:
less postoperative adjustments
- 3 inner compartment
integrated in band:
precise and concentric filling



first closing position 10.5 cm (\varnothing 31 mm)
second closing position 9.75 cm (\varnothing 26 mm)

precise and concentric filling



MiniMIZER EXTRA



- 4 pre-shaped gastric band
for easy positioning
- 5 small port (\varnothing 13 mm)
for sub-xyphoidal placement
- 6 medium port (\varnothing 18 mm)
for sub-costal or rectal sheath placement



MiniMIZER EXTRA



with loops



 **BARIATRIC
SOLUTIONS**

Bariatric Solutions GmbH
Kaltenbacherstr. 24
CH - 8260 Stein am Rhein

Phone: +41 - 52 741 - 65 00
Fax: +41 - 52 741 - 65 02
info@bariatric-solutions.com
www.bariatric-solutions.com